

Multi-Calibrators™

I125 & Co57 Sources

- Superior accuracy
- Most economical
- Best matching
- Regulatory documentation
- Highest quality
- Most experience
- Since 1985



See reverse side for product specifications

Bind-It™ Decontaminant

Radioactive Iodine Remover for All Surfaces

- Removes all radioactive iodine (I123, I125, I126, I131)
- Extremely effective
- Safe for Gamma Counter & Camera Detectors
- Super-strong ionic bond
- Very affordable



See reverse side for product specifications

Bind-It™ Handsoap

Radioactive Iodine Remover for Hands

- Special emollients leave skin soft
- Anti-Bacterial
- Removes radioactive iodine from skin
- Extremely effective
- Gentle to skin
- Very affordable



See reverse side for product specifications

Multi-Calibrators™

The **Gold Standard** by which all other I125 & Co57 calibrators are judged. Since 1985, LTI has been supplying the industries' premiere gamma counter calibrators. Due to our high volume and attention to detail, our calibrators are the highest quality available. This high volume also allows us to be the most affordable on the market.

* As the world's leader in gamma counter technology, we know how best to calibrate them.

LTI Multi-Calibrators™ feature:

Multi-Log™: Complete records for Regulatory compliance.

Life Span: I125 last for 6 months & Co57 for 18 months

Superior matching: $\leq 0.5\%$

High DPM: We ship all calibrators with a $DPM \geq 150,000$

Bind-It™ Decontaminant

All forms of radioactive iodine are highly volatile. This is precisely why they make such excellent diagnostic labels. However, this volatility if not properly dealt with, can cause contamination of: gamma counters, camera's, uptake probes, pipettes, lab benches, furniture, lab coats etc. Bind-It™ decontaminant has been specially formulated to remove contamination from all surfaces. It is safe and effective for use on delicate gamma counter/camera detectors. It is noncorrosive and safe to skin. Special anti-bacterial properties further aid in keeping your laboratory safe and clean.

Because Iodine is uniquely negative in its electrical charge, other radioactive removers rarely work; they may actually push iodine away, causing the contamination to spread. Only Bind-It™ is specially formulated to remove radioactive Iodine.

Bind-It™ Handsoap

Bind-It™ Handsoap uses the same active ingredients as Bind-It™, but in a gentle handsoap. Added emollients leave your hands soft and moisturized, while ridding them of unwanted radioactive Iodine and bacteria. Gentle enough for routine use, this handsoap is economical enough to replace your regular laboratory handsoap while giving an extra measure of protection from contamination and infection.

*LTI manufactures the extensive Genesys™ line of single & multi-well gamma counters which are sold worldwide.

Specifications:

Multi-Log™

Complete, detailed regulatory records.

Life Span

I125 lasts 6 months, Co57 for 18 months

Best matching

$\leq 0.5\%$

Shorter Count Time

Initial DPM $\geq 150,000$

Bind-It™ Concentrate

16 oz. concentrated cleaner

Dilute with H₂O at 20:1 to use. **Economical**

Bind-It™ Decontaminant Spray

32 oz. high volume spray that helps to lift iodine from surfaces. **Ready to Use**

Bind-It™ Decontaminant Wipes

Ready to use Moistened Towelets

in an easy to use dispenser. **Convenient**

Bind-It™ Handsoap

(12) 32 oz. pump bottles

(1 or 4) U.S. Gallon pump bottle



43W900 Route 64
Maple Park, IL 60151
800.542.1123
630.365.1000
630.365.9687 Fax
www.labtechinc.com